

Hands-On Lab: Import and Export project as a Solution

Scenario


You received an APMS solution to train how to import and export it.

Lab Objectives

- Import the file TrainingProjectSolution.apms
- Export the recently imported solution as a private solution

Steps


1. Make sure you have downloaded the file **TrainingProjectSolution.apms** from the learning portal
2. Go to the Administration area of ARIS Process Mining



3. In the administration section go to **Install solution** and either drag and drop the file or click on Select file and choose the provided .apms file
 - a. Check the information of this solution and click on install
 - b. In the Configure step you may choose which solutions assets you would like to import. Press Next
 - c. Now you may choose to create a new Project or overwrite an existing project. Create a new project with this information

- i) Name: Training Solutions
- ii) Abbr: TS
- iii) Color: Your choice

4. Create a new Data Set with the same name of the project. Press next to install the solution.
5. You may check the project (Go to Project). This project as stated in the description contains the solution of other hands on labs. Feel free to consult it when needed.
6. To export you project as a solution go back to the Projects section.
7. Click on the three dots near the Training Solution and choose **Export as Solution**



8. Select your solution type Private and if you want you may also change the name of the solution
9. In the configure part select all assets.
10. And just click on the download button in the end to download our solution